

Peaks table

Run

Bank

UB index

Pred\Found

Sat\Main

Run number	DetID	h	k	l	Wavelength	Energy	TOF
0	201	0	0	0	1.99779684	20.496	
0	101	0	0	0	2.28319550	15.692	
0	102	0	0	0	2.28319174	15.692	
0	203	0	0	0	1.99779428	20.496	

Table can be sorted by any column including by I over sigma (default)

4 out of 123 peaks

Select AllSelect All Filtered

LoadFind PeaksFind UBChoose CellCalibrateSample InfoConvert to QIntegrateExportVisualise

Run Number(s)

SXD33304

Browse

V run

Browse

Empty instrument

Browse

Goniometer

X	Y	Z	Handedness	Angle\Log name
0	1	0	+1	CCR_Pos
0	0	1	+1	45
0	1	0	-1	EWALD_Pos

Load

Calibration File

SXD_NaCl_cycle19_3.xml

Browse

Apply

Run	V Normalised	Calibrated
033304	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
033305	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
033306	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Remove SelectedRemove All

Log output

Loading geometry cache from C:\Users\xxxxx\AppData\Local\mantid\instrument\geometryCache\SXDadd3e2457fedeceb8814ed6d2b4aaf0c3ee2\Parsing from XML file: C:\mantid\instrument\SXD_Load successful, Duration 5.47 seconds

CopySelect AllClear AllScrollbar limitLog Level>

UI would need to remember the last used V run to avoid user having to respecify on each load

NormaliseByCurrent on Load as well

Goniometer log names remembers last ones used

calibration .xml file

Removing a run from the runs table deletes any associated peaks from the peaks table also

Each tab has a log output window showing output from processes initiated on that tab - so no need to go to master workbench message log

[illegible]

Args should be populated with the default value for all mandatory arguments

Populate peak markers on the instrument view (Visualise tab)
No slice viewer at this point

Peaks table gets populated

[illegible]

[illegible]

[illegible]

Load Find Peaks Find UB Choose Cell Calibrate Sample Info Convert to Q Integrate Export Visualise

Shape

Sphere

Mesh File

Browse

Dimensions

0.000564m

Chemical Formula

C3-H7-N-O2

Number Density

0.0093591

Apply

Log output:

4 out of 123 peaks

Select All

Select All Filtered

Run SetSample on all workspaces (all matrix and peaks workspaces)

Shape drop down has:

- Sphere
- Cylinder
- Mesh

If Mesh selected the Mesh File Browse box is enabled

Dimensions box will be shape specific eg radius only for sphere, radius, height, axis for cylinder etc

For mesh and maybe other shapes, need to supply initial orientation. May imply extra fields

Need to store the V sample details somewhere but won't vary for same instrument so put in settings somewhere?

[illegible]

Users may select runs by multiselect run dropdown or specifying UB index
Only permit one selected UB index if output dimensions=HKL
Args should be populated with the default value for all mandatory arguments
Use algorithm: ConvertToDiffractionMDWorkspace
Slice viewer and 3D isosurface view activate (Visualise tab)

Output Dimensions:

- Q (sample)
- Q (lab)
- HKL

Possibly in future add k_f to possible output dimensions

May be more generic to allow bank selection via the instrument tree in instrument viewer rather than a bank drop down. Bank definition varies (level/type) across instruments

[illegible]

Available algorithms:

- Integrate in Detector space (IntegratePeaksSkew + perhaps others in future)
- Integrate in QLab (IntegratePeaksMD)
- Integrate in HKL (IntegratePeaksMDHKL)

Args should be populated with the default value for all mandatory arguments

The set of peaks/runs that are to be integrated is driven by the peaks table on the left hand side. If the user selects peaks from more than one UB then the IntegrateHKL algorithm should be disabled

IntegratePeaksMD should be disabled if there isn't an MD workspace for every run that's been selected

Predicted\Found filter should never allow both to avoid integrating huge number of peaks

[illegible]

